


Fact Sheet

UniMedCar10

telelift
SIMPLIFYING INTRALOGISTICS



Safe and hygienic transports in hospitals

The UniMedCar10 stands for a safe and hygienic transport in hospitals.

It transports goods, up to 10 kg transport weight, quickly and reliably between operating theatres, pharmacies, laboratories, storage areas and wards in hospitals.

All installed components in the UniMedCar10 stand for a long service life and maximum reliability.

Antimicrobial surface coating

Particularly hygienic transport is ensured by an antimicrobial surface coating, which prevents the growth and spread of many bacteria and germs on the surface of the trolley.

Takes over the transport logistics for:

- Samples (Single- and bulk transport)
- Medication (Single- and bulk transport)
- IV-bags and blood plasma
- Infusions
- Vaccines
- Sterile goods
- Radioactive goods with safe and secure shielding
- Documents and media

Features UniMedCar10:

- Up to 10 kg transport weight
- Max. speed: 0,6 m/s horizontal
- Antimicrobial surface coating
- Optional with internal swivel box
- Electronic lid closure
- Safety buffer as contact protection

Fact Sheet

UniMedCar10

telelift
SIMPLIFYING INTRALOGISTICS

Transportgewicht und -abmessungen:



Ca. 38 Litres
Max. Transport volume



503 mm x 185 mm x 415 mm
Internal dimensions (l x w x h)



max. 10 kg
Transport weight



5° to 40° C
Operation temperature

